



# The cutting edge solution to detect donor-derived cell-free DNA for transplant surveillance

**This transplant focused innovation in organ transplant surveillance can drive better outcomes for your patients**

AlloSeq cfDNA is a kit\* or service\*\* based solution that enables measurement of dd-cfDNA through a blood test




\* AlloSeq cfDNA kit is CE-IVD in EU and the UK. In the rest of the world, AlloSeq cfDNA kit is for Research Use Only. Not to be used for diagnostic procedures.

\*\* AlloSeq cfDNA Service is for Research Use Only. Not to be used for diagnostic procedures.


AlloSeq cfDNA is only available outside of the United States.



# Cell-free DNA: a clear biomarker for organ injury




## AlloSeq cfDNA, innovation in action




# What is AlloSeq cfDNA?

- + A minimally invasive blood test that measures dd-cfDNA
- + AlloSeq cfDNA kit\* is an NGS based test to run in your own lab and does not require prior genotyping
  - A simple workflow with sample to report solution within 24 hrs
- + AlloSeq cfDNA Service\*\* is an alternate way to send your samples into CareDx lab to validate or use AlloSeq cfDNA without needing NGS equipment
  - Simple sample collection process with an easy-to-interpret report delivered directly to you



AlloSeq cfDNA can be added to these current protocols with ease:



\* AlloSeq cfDNA kit is CE-IVD in EU and the UK. In the rest of the world, AlloSeq cfDNA kit is for Research Use Only. Not to be used for diagnostic procedures.

\*\* AlloSeq cfDNA Service is for Research Use Only. Not to be used for diagnostic procedures.

AlloSeq cfDNA is only available outside of the United States.

# AlloSeq cfDNA was developed after the successful clinical validation of AlloSure


AlloSure is the dd-cfDNA surveillance service that has been clinically and analytically validated for identifying kidney, heart and lung injury.<sup>1-6</sup>

AlloSure is the most published dd-cfDNA test in transplantation with over 35+ peer reviewed journal publications.


## Clinical Decision Support with AlloSure

AlloSure at varying levels provides meaningful, actionable information to your existing clinical paradigm


Interpretation of AlloSure dd-cfDNA in Kidney Transplantation<sup>1-4</sup>



Interpretation of AlloSure dd-cfDNA in Heart Transplantation<sup>5\*</sup>




Interpretation of AlloSure dd-cfDNA in Lung Transplantation<sup>6\*</sup>



# AlloSure Kidney: Data Driven Insights

AlloSure is highly sensitive in distinguishing ABMR from no ABMR



**AlloSure > 2.9%**


**89% PPV** for ABMR  
in DSA+ Patients<sup>1</sup>

Sensitivity: 50%  
Specificity: 94%  
Prevalence: 40%

} At 2.9% dd-cfDNA

**2.9% is the median from DART patients with ABMR in DSA+ Patients**

AlloSure detects active rejection




**AlloSure > 1.0%**

**threshold**

for Active Rejection<sup>\*2</sup>


AlloSure differentiates biopsy confirmed TCMR1A / Borderline rejection patients



**AlloSure > 0.5%**  
**for eGFR Decline**

in TCMR1A/Borderline Patients<sup>3</sup>

Reassurance with a low AlloSure score



**AlloSure < 0.21%**

**has 95% NPV**

for Active Rejection<sup>\*4</sup>

Sensitivity: 85%  
Specificity: 33%  
Prevalence: 10%

} At 0.21% dd-cfDNA

**0.21% is the median from DART healthy stable recipients**

\*Active Rejection = acute/active ABMR; chronic, active ABMR; and TCMR IA and greater

# dd-cfDNA surveillance testing can be added to previously established protocols

## Example - For a Kidney Transplant:

### Months 1, 2 post-transplant

- + Baseline to be established during this time
- + dd-cfDNA is associated with changes immediately post-transplant (ischemia reperfusion injury, DGF or medication dose changes) and may be a potential surrogate marker for these changes
- + Results may contribute to treatment decisions

### Months 3, 4 post-transplant

- + dd-cfDNA changes can be used as an associated surrogate marker to monitor the waning of induction immunosuppressants
- + Continue to monitor changes from baseline
- + Maintain continuity in surveillance for patients who transition to a general nephrologist

Biomarker	Condition Tested	Time post-transplant	
		1 week	1 Month
Creatinine	Indirect graft function	Daily	2-3 per Week
BK Virus	Viral infection	Monthly	
DSA	Donor Specific Antibody formation	Weekly	Monthly
<b>dd-cfDNA</b>	<b>Active allograft injury</b>	<b>Monthly</b>	



### Months 6, 9, 12 post-transplant

- + Monitor for acute rejection
- + Test along with DSA to improve the positive predictive value for antibody mediated rejection<sup>1</sup>

### Year 2 Onward: quarterly testing through life of transplant

- + Monitor for changes in dd-cfDNA results that may occur before the onset of symptoms
- + Rejection due to ABMR increases over time and monitoring may provide an early warning of rejection<sup>1</sup>
- + A significant increase in dd-cfDNA may provide early insight into medication non-adherence

2-3 Months	4-6 Months	7-12 Months	12+ Months
Weekly	Every 2 Weeks	Monthly	Every 2-3 Months
	Every 3 Months		-
-	Every 6 Months		
	Every 2 Months	Every 3 Months	

## ORDERING INFORMATION

Product	No. RXNs	Product Code
AlloSeq cfDNA - IVD	24	ASCF.1(24)-IVD
AlloSeq cfDNA - RUO	24	ASCF.1(24)
AlloSeq cfDNA Software	N/A	ASCFS1.0
AlloSeq cfDNA Service	N/A	N/A

AlloSeq cfDNA kit is CE-IVD in EU and the UK. In the rest of the world, AlloSeq cfDNA kit is for Research Use Only. Not to be used for diagnostic procedures.

AlloSeq cfDNA Service is Research Use Only. Not to be used for diagnostic procedures.

AlloSeq cfDNA is only available outside of the United States.

## References

1. Jordan SC et al. *Transplant Direct*. 2018; 4:e379
2. Bloom RD et al. *J Am Soc Nephrol*. 2017; 28:2221–2232
3. Stites E, et al. *Am J Transplant*. 2020; 00:1–8
4. Bromberg JS et al. *J Appl Lab Med*. 2017; 2:309–321
5. Khush KK, et al. *Am J Transplant*. 2019 Oct;19(10):2889-2899; D-OAR is a sub-study of the Outcomes AlloMap Registry (OAR)
6. Keller, M.B et al; *The Journal of Heart and Lung Transplantation*, 40(4), S148



### Want to order AlloSeq cfDNA surveillance tests?

- **AlloSeq cfDNA Services:** Send in samples to CareDx central lab in Stockholm, Sweden, Email: [servicetesting@caredx.com](mailto:servicetesting@ caredx.com) or call +46-8-50893900

- **AlloSeq cfDNA Kit:** To learn more about how to bring AlloSeq cfDNA to your own lab, contact your CareDx representative or reach out to us:

Americas [orders-US@caredx.com](mailto:orders-US@caredx.com)  
EMEA [orders-se@caredx.com](mailto:orders-se@caredx.com)  
APAC [orders-aus@caredx.com](mailto:orders-aus@caredx.com)

AlloSeq cfDNA kit is CE-IVD in EU and the UK. In the rest of the world, AlloSeq cfDNA it is for Research Use Only. Not to be used for diagnostic procedures.

AlloSeq cfDNA Service is Research Use Only. Not to be used for diagnostic procedures.

AlloSeq cfDNA is only available outside of the United States.



The cutting edge solution to detect donor-derived cell-free DNA for transplant surveillance

[www.caredx.com/alloseq-cfdna](http://www.caredx.com/alloseq-cfdna)