



SPECI-GARD®



The bag SPECI-GARD® is used as internal hospital and external secondary packaging for the transport of infectious substances of Category B (**UN3373**), compliant to International and European standard and guideline for the transport of infectious substances. (**ADR2013/P650 - WHO/HSE/GCR/2012.12**).

TECHNICAL DATA

Specification

Material

Thickness

Size (±5 mm)

Closing system
Opening system

Structure of the bag

Printing

In option

Packing
Storage
Expiry date

Information

Origin

DESCRIPTION

The SPECI-GARD® is a **liquid-tight** bag in polyethylene MICRO-TEX® with adhesive closure system. Qualified liquid-tight by the LNE and the B.V.T according to the instructions of packaging **P650 of the class 6.2 of the ADR**.

Polyethylene low density - 100% Virgin, 100 % recyclable
- Low environment impact

40 µm (±5 %)

165 x 254 mm - ref. : 95600BIO
185 x 254 mm - ref. : 95700BIO
254 x 254 mm - ref. : 95610BIO
254 x 335 mm - ref. : 951013BIO

Secure, **permanent adhesive closure**

Simple and fast opening thanks to a unique and innovative precut (not requiring sharp or cutting objects).

Liquid-tight pouch dedicated for the specimen sample.
Additional pouch for the inquiry/requisition form.

Logo Biohazard (infectious risks) - Operating Instructions - Number of lot.
Ink used for the printing of bags is guaranteed unleaded and without heavy metals.

- Colors red, green, blue, yellow, orange, etc.
- System of head of bundle for wall suspension of bags
- Customized printing (logo, name, address, color)
- Pad absorbent inserted or separately

2 000 bags per carton.
Ambient temperature.
Better use before 3 years.

Delivery 48 to 72 hours for standard bags or 6 to 8 weeks for customized bags.
Contact us for any customization (printing and size).
Test report of liquid-tightness by external laboratories on request.

France

