

New From Aferrix

feROS™ **eLPI** Ref. TSE804

An Enhanced Labile Plasma Iron (**eLPI**) test for the simultaneous detection of two different forms of Non Transferring Bound Iron (NTBI): overt and cryptic.

Overt LPI - is the redox active and chelatable iron fraction of NTBI detected by the feROS™ LPI test.

Cryptic LPI – is the redox inactive fraction of NTBI which is revealed by the feROS™ test in the presence of a mobilizer.

feROS™ eLPI – for comprehensive detection of NTBI.

feROS™ **eLPI** is available in kit format: Ref. TSE904

For more information, please contact us at:
info@aferrix.com



Aferrix Ltd.

19A Habarzel St, Ramat Hachayal, Tel Aviv 69710, Israel
Tel: 972.73.2126155 - Fax: 972.73.2126144
info@aferrix.com
www.aferrix.com